



evercare® Basic Isolation Gown with thumb loops

REF: 270009, 270009-1, 270009-2, 270009-3, 270010, 270010-1, 270010-2

evercare® Basic Isolation Gown with thumb loops is an impervious isolation gown with thumb loops made of fluid resistant polyethylene to protect clothing from spills and splashes during medical and cleaning procedures.

Key Features:

- ▶ Open back design
- ▶ Thumb loops sleeves to keep sleeves from sliding up when donning gloves
- ▶ Quick to take on and off by tearing the neck section
- ▶ Individually folded and packed in unit pouch



REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
270009	evercare® Basic Isolation Gown with thumb loops, light blue	104 x 120 cm	25/100
270009-1	evercare® Basic Isolation Gown with thumb loops, light blue	98 x 135 cm	25/100
270009-2	evercare® Basic Isolation Gown with thumb loops, light blue	93 x 112 cm	25/100
270009-3	evercare® Basic Isolation Gown with thumb loops, light blue	94 x 114 cm	25/100
270010	evercare® Basic Isolation Gown with thumb loops, light blue	110 x 140cm	25/100
270010-1	evercare® Basic Isolation Gown with thumb loops, light blue	112 x 117 cm	25/100
270010-2	evercare® Basic Isolation Gown with thumb loops, light blue	114 x 135 cm	25/100

Sales Unit: Transport Carton

Technical Summary

Regulations

The product complies with legal requirements for medical devices



REACH regulation:	The product(s) covered by this data sheet does(do) not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Compliance with legal requirements:	Medical Devices Directive 93/42/EEC, Class I
Compliance with product standards:	N/A
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN 1041

Material

Product Part	Raw Material & Characteristic
Main material	Cast polyethylene film (CPE), light blue, embossed finish.

Warnings and Recommendations

Warnings:	The product is single use. If reused performance of product may deteriorate, cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 Years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	PE wallet pouch	To be confirmed
Dispenser Pack	PE bag	To be confirmed

Transport Carton	Carton, TCF quality	400 x 200 x 350
Packing Tape	Transparent, PVC-free	N/A

Barcodes

REF	Single pack	Dispenser	Transport Carton
270009	N/A	 6 438129 324525	 6 438129 524529
270009-1	N/A	 6 438129 324532	 6 438129 524536
270009-2	N/A	 6 438129 324556	 6 438129 524550
270009-3	N/A	 6 438129 324570	 6 438129 524574
270010	N/A	 6 438129 324495	 6 438129 524499
270010-1	N/A	 6 438129 324549	 6 438129 524543
270010-2	N/A	 6 438129 324556	 6 438129 524550

ONEMED

Manufactured in China for
OneMed Group Oy
 Address: Metsälantie 20
 Postal Address: P.O. Box 10
 FI-00321 Helsinki, Finland
 Phone +358 20 786 6800
www.onemed.com
 Domicile Helsinki
 Business ID 2039640-1